

SAP ABAP table SLS_XET {PAW - Tests per Event}

SLS_XET
<u>MANDT</u> [CLNT (3)]
PLVAR [CHAR (2)]
<u>EVEID</u> [NUMC (8)]
<u>NO_TEM</u> [CHAR (1)]
<u>TESTID</u> [NUMC (8)]
<u>TTYID</u> [NUMC (8)]

SLS_EVENT
<u>MANDT</u> [CLNT (3)]
PLVAR [CHAR (2)]
<u>EVID</u> [NUMC (8)]
<u>NO_TEM</u> [CHAR (1)]

SLS_TEST
<u>MANDT</u> [CLNT (3)]
<u>TESTID</u> [NUMC (8)]

T000
<u>MANDT</u> [CLNT (3)]

